

This would be a screen-by-screen menu just like on the Somnormal app

1. First-time users should register first.

SCREENSHOT OF THE REGISTRATION SCREEN

2. Select **Pairing** from the drop-down menu and press **Search for Cervella** to enter the pairing menu.

SCREENSHOT OF THE HOME MENU WITH THE PAIRING MENU HIGHLIGHTED

3. Press **“Search Again.”** If your Cervella is nearby and charged, it should appear in the list below. Select “Cervella-XXXX” and the Cervella LED should illuminate green. You have successfully paired the Cervella to your smart device.

SCREENSHOT OF THE PAIRING MENU WITH THE CERVELLA IDENTIFIED

4. Apply thin layer of conductive gel or (if not available) water to the gold-colored portion of the ear cushions and at the point of skin contact behind the ear area (mastoid) as circled below:



PHOTO OF APPLYING GEL TO EARMUFF CUSHION

5. Connect the Cervella device to the headset using the white connection cable.



PHOTO OF CONNECTING CERVELLA TO HEADSET

6. Place the headset on your head ensuring that the ear cushions are making tight contact with your skin, especially in the areas where the conductive gel was applied.



PHOTO OF HEADSET ON HEAD WITH GOOD SKIN CONTACT